

SILVRSTAT® ANTIBACTERIAL WOUND DRESSING GEL

Manufactured for ABL Medical, LLC American Fork, UT 84003 (801) 763-8000

PRODUCT DESCRIPTION

SILVRSTAT® Antibacterial Wound Dressing Gel is a water-based gel wound dressing for use in the management of first and second degree burns, wounds such as stasis ulcers, pressure ulcers, diabetic ulcers, lacerations, abrasions, skin tears, surgical incision sites, device insertion site wounds, graft sites and donor sites.

SILVRSTAT® Antibacterial Wound Dressing Gel contains silver that in laboratory tests has been shown to inhibit the growth of microorganisms such as: *Staphylococcus aureus*, *Pseudomonas aeruginosa*, *Escherichia coli*, antibiotic resistant strains, *MRSA* and *VRE*, as well as fungi such as *Candida albicans*, and *Aspergillus niger*.

Indications for Use

SILVRSTAT® Antibacterial Wound Dressing Gel is indicated for the management of:

- First and second degree burns
- Wounds such as: •stasis ulcers •pressure ulcers •diabetic ulcers
- Lacerations •Abrasions •Skin tears
- Surgical incision sites
- Device insertion site wounds •Graft sites
- Donor sites

Contents: 32 PPM Proprietary Silver (purified water, nano-silver at .01 micron), Propylene Glycol, Triethanolamine, Carbomer

Directions for Use:

SILVRSTAT® Antibacterial Wound Dressing Gel should be applied directly to the wound and then covered with a conventional dressing. SILVRSTAT® Antibacterial Wound Dressing Gel should normally be changed every 24 to 72 hours, but the frequency will depend on the type and condition of the wound.

Precautions: Intended for external use only. **If your condition persists or gets worse, contact your healthcare provider.**

Adverse Reactions: There are no known or reasonably suspected adverse reactions associated with the use of SILVRSTAT® Antibacterial Wound Dressing Gel.

Storage: Store at room temperature.

Caution: Federal law restricts this device to sale by or on the order of a licensed healthcare practitioner.

